

## INSTANT-VIEW<sup>®</sup> FECAL OCCULT BLOOD TESTS



- Sensitivity: 50 ng hHb/ml buffer or 50 µg hHb/g feces
- Specificity: Over 98% specific to human hemoglobin with no dietary restrictions
- Results in 5 to 10 minutes
- FDA 510(k) cleared
- CE Certified for sale in European Union
- CLIA Waived

*Rapid, one-step immunochemical test intended for the qualitative detection of Fecal Occult Blood by laboratories or physicians offices.*

## INSTANT-VIEW<sup>®</sup> PROSTATE SPECIFIC ANTIGEN (PSA) TESTS

- Sensitivity: 4 ng/ml and/or 10 ng/ml
- Results in 4 to 7 minutes
- Available in serum only or serum/whole blood combination
- For export only, not for sale in the U.S.



*Initial qualitative screening test designed to detect prostate specific antigen (PSA) concentration in serum or whole blood.*

00-0104 REV A

**ALFA SCIENTIFIC DESIGNS, INC.**

FDA Registered • ISO Certified

TEL: 858.513.3888 / Toll Free (U.S. only): 877.204.5071 / FAX: 858.513.8388

WWW.ALFASCIENTIFIC.COM

