



ALFA NEWSLETTER

Making Drug Testing Really Simple

Making its debut this spring at Alfa Scientific Designs, Inc. is the One-Step Simple DOA Cup.

Now available in both generic and Instant-View® labeling, the Simple DOA Cup combines both a sleek one step design and the CLIA Waived strips that Alfa manufactures. While not all potential configurations are waived, this does add to the menu of cup selections available from Alfa.

Alfa is always striving to bring you the latest and greatest products on the market. The Simple DOA Cup makes an excellent addition to our portfolio, joining Alfa's patented Split-Screen cup design.

Initially we will be

offering the following popular configurations in one or both of the cup formats:

5 Panels

- AMP, COC, MET, MOR, THC (CLIA Waived)

- AMP, COC, MOR, PCP, THC (CLIA Waived)

- COC, MET, MOR300, PCP, THC

6 Panel

- AMP, BZD, COC, MET, MOR, THC (CLIA Waived)

12 Panel

- AMP, BAR, BZO, BUP, COC, MET, MDMA, MOR, MTD, OXY, PCP, THC

All of the above listed combinations will include adulteration testing with parameters of Nitrite, Creatinine, pH, Specific Gravity, Glutaraldehyde, Oxidant.

For more information about the Simple or Split Screen Cups please contact your Alfa representative today.



Your Source For Reliable, Accurate, Rapid Diagnostic Testing Products for Fertility, Drugs of Abuse, & Infectious Disease
MADE IN THE U.S.A.

Check out www.alfascientific.com



alfa

SCIENTIFIC DESIGNS, INC

13200 Gregg Street, Poway, CA, 92064

PH: 858-513-3888 FX: 858-513-8388

Volume 1, Number 2 - Spring 2008 (Issue 2)

Editor: Alison Adduono aadduono@alfascientific.com

3-in-1 Cardiac Test Coming Soon

By RICK GUNNILL

Alfa has filed for US FDA 510(k) clearance to market its new Cardiac 3-in-1 Whole Blood/Serum Test (cTnI, CKMB, MYO) test. The product is intended for professional use as an aid to the rapid diagnosis of acute myocardial infarction (AMI).

Alfa's Cardiac 3-in-1 test, which is available in cassette and dip-strip formats, can use either whole blood or serum samples. Each device has a test strip with separate test lines for each of three markers: cardiac Troponin I, creatine kinase MB isoenzyme (CKMB) and myoglobin. Cut-off concentrations for these markers in Alfa's test (1.5 ng/ml

Free cTnI, 5 ng/ml CKMB and 100 ng/ml myoglobin) indicate when their levels have risen above normal baseline levels. In humans, concentrations of the markers in blood increase at different rates after an AMI event to reach their respective maxima at different times.

Thus, in comparison with a single marker, testing the 3 markers can be more informative to a physician evaluating symptoms to make a diagnosis regarding AMI.

Stay tuned for more information on this new product in coming months.



Alfa's cardiac marker test devices.

Alfa Scientific is attending the following tradeshow this year.
Stop by our booth... we'd love to see you.



GSA Expo
April 22-24
Anaheim, CA
Booth # 119
<http://www.expo.gsa.gov>



Datia
April 30-May 1
Dallas, TX
Booth # 404
<http://datia.org/>



AACC
July 27-31
Washington, DC
Booth # 156 & 158
<http://www.aacc.org>



Medica
Nov. 19-22
Germany
Booth # TBA
<http://www.mdna.com>



Corporate Web Site Gets Facelift

If you haven't already noticed, Alfa Scientific's corporate web site has undergone a dramatic transformation in recent months. Not only are some new product offerings being featured but it has a completely new look and feel.

One of the goals of our revamp is to make it a vital resource for our customers - a source of product information as well as helpful tools such as drug detection ranges, certification quizzes, newslet-

ters, etc. In coming weeks you will see these items added to our website so please check it out from time to time. We're also adding an On Line Store for the convenience of our customers and to make over-the-counter drug testing supplies available.

On Line Store

Coming soon is Alfa Scientific Designs' 24/7 On Line Product Store where we will feature our most popular

tests for at home use. You can order as many tests as you like all from the comfort and privacy of your home. All products are shipped out in a discreet packaging within 24 hours of receiving your order Monday thru Thursday.

Please visit our corporate web site at www.alfascientificdesigns.com to learn more about the company and product line.

State Laws

Did you know that many states have their own laws and regulations that impact when, where and how employers can implement their drug-free workplace programs?

Check out this website for more detailed information:
<http://www.dol.gov/asp/programs/drugs/said/StateLaws.asp>

ALFA
SCIENTIFIC DESIGNS, INC.

Quality Made in America

Home | Alfa Products | Product Spotlight | Become a Distributor | News & Events | Product Information | FAQs | Buy Products Online

Quality, easy-to-use products.

ALFA Scientific Designs, Inc. develops and manufactures high quality rapid in-vitro diagnostic devices to serve the point-of-care marketplace including hospitals, physicians, criminal justice, and occupational health centers. We focus on rapid immunoassays including those for the detection of drugs of abuse, fertility hormones, infectious diseases, cardiac and cancer markers. ALFA is committed to satisfying customer's needs by ensuring the highest level of performance of every device we manufacture.

Quality. Service. Commitment.

GSA Advantage!
View Products
Contract # GS-07F0428T
GSA Members order contact:
mlathe@alfascientific.com
(858) 513-3888

AACC 2008 AACC Annual Meeting
July 27 - 31, 2008
Washington, DC

MEDICA 2008
November 19 - 22, 2008
Dusseldorf, Germany

DATIA DATIA Annual Conference 2008
April 30 - May 1, 2008
Dallas, Texas

2008 GSA EXPO 2008 GSA Expo
April 22 - 24, 2008
Anaheim, CA

INSTANT-VIEW®
Easy to Read
Split Screen Cup
New!

Alfa Instant-View® Drug of Abuse Test

Terms Privacy Customer Service Technical Support Company Profile Contact 13200 Gregg Street, Poway, CA 92064 USA • Call 877-204-5071

Pictured above is a snapshot of the new home page for Alfa's corporate website.

