

INSTANT-VIEW[®] Simple, Fast Testing To Help Screen For Cardiac Events



- Simple One-Step Test
- Over 94.2% Specific
- Sensitivity:
 - 1.5 ng/mL free cTnI or 5 ng/mL Complex cTnI
 - 5 ng/mL CKMB
 - 100 ng/mL Myoglobin
- Results in 15 Minutes
- Convenient Room Temperature Storage
- Made in the USA

Product Description	Kit Size	FDA Cleared	CE Mark
Troponin I (TNI) Serum Cassette	25	✓	✓
Troponin I (TNI) Whole Blood/Serum Cassette	25	✓	✓
Troponin I, CKMB, Myoglobin (3-in-1) Whole Blood/Serum Cassette	25		✓

To order or for further information please contact one of our representatives today:



13200 Gregg Street, Poway, CA 92064

PH: 858-513-3888 • TOLL FREE: 877-204-5071 • FAX: 858-513-8388

WEB: www.alfascientific.com • EMAIL: info@alfascientific.com