

Simple, Fast Testing For Detecting Drugs of Abuse



- FDA 510(k) for Professional Use
- CE Marked for Sale in Europe
- OTC Cleared and CLIA Waived



INSTANT-VIEW[®] DOA Tests

One-Step Drugs of Abuse Immunoassays

For the detection of illegal or prescription drugs in human urine to determine illicit use or pain management adherence.



WHERE VALUE MEETS QUALITY IN DIAGNOSTICS™

INSTANT-VIEW[®] Simple, Accurate DOA Testing

92.2%-99.8% Accurate • Results in 4-7 Minutes • Test up to 12 Drugs at Once

Made in the USA • Adulteration Controls Available • Dip Strip and Device Formats

DRUG ANALYTE CHOICES					
Product	Description	Cut Off	FDA Cleared	Waived	CE Mark
AMP	Amphetamine	1000 ng/mL	✓	✓	✓
AMP300	Amphetamine	300 ng/mL	✓		✓
BAR	Barbiturate	200 ng/mL	✓	✓	✓
BZD	Benzodiazepine	300 ng/mL	✓	✓	✓
BZD200	Benzodiazepine	200 ng/mL			
BUP	Buprenorphine	10 ng/mL	✓		✓
COC	Cocaine	300 ng/mL	✓	✓	✓
COC150	Cocaine	150 ng/mL	✓		✓
EDDP	Methadone Metabolite	100 ng/mL			
KET	Ketamine	1000 ng/mL			
MET	Methamphetamine	1000 ng/mL	✓	✓	✓
MET500	Methamphetamine	500 ng/mL	✓		✓
MET300	Methamphetamine	300 ng/mL	✓		✓
MOR/OPI	Morphine/Opiate	2000 ng/mL	✓	✓	✓
MOR/OPI300	Morphine/Opiate	300 ng/mL	✓		✓
MTD	Methadone	300 ng/mL	✓	✓	✓
OXY	Oxycodone	300 ng/mL	✓		
OXY100	Oxycodone	100 ng/mL	✓		✓
PCP	Phencyclidine	25 ng/mL	✓	✓	✓
PPX	Propoxyphene	300 ng/mL	✓		✓
TCA	Tricyclic Antidepressant	1000 ng/mL	✓	✓	✓
THC	Marijuana	50 ng/mL	✓	✓	✓
TRA	Tramadol	200 ng/mL			✓
XTC	Ecstasy/MDMA	500 ng/mL	✓	✓	✓

To order or for further information please contact one of our representatives today:



13200 Gregg Street, Poway, CA 92064

PH: 858-513-3888 • TOLL FREE: 877-204-5071 • FAX: 858-513-8388

WEB: www.alfascientific.com • EMAIL: info@alfascientific.com

00-0117 F

WHERE VALUE MEETS QUALITY IN DIAGNOSTICS™