

Simple, Accurate Testing To Help Screen For Menopause



- FDA 510(k) Exempt for OTC or Professional Use
- CE Marked for Sale in Europe



INSTANT-VIEW[®] FSH Test

One-Step Menopause Immunoassay

For the detection of follicle stimulating hormone (FSH) in human urine to aid in the diagnosis of menopause.



WHERE VALUE MEETS QUALITY IN DIAGNOSTICS™

INSTANT-VIEW[®] Simple, Accurate Testing To Help Screen For The Onset Of Menopause



- Simple One-Step Test
- Over 99% Accurate
- Fast Results in 5 Minutes
- Convenient Room Temperature Storage
- Diagnose Onset of Menopause
- Multiple Test Formats
- Made in the USA

Product Description	Kit Size	Cut-Off Level	FDA Cleared	Waived	CE Mark
FSH Urine Cassette	25	40 mIU/mL	Exempt		
FSH Urine Dip Strip	50	40 mIU/mL	Exempt		
FSH Urine Midstream (OTC Use)	1	40 mIU/mL	Exempt		

To order or for further information please contact one of our representatives today:



13200 Gregg Street, Poway, CA 92064

PH: 858-513-3888 • TOLL FREE: 877-204-5071 • FAX: 858-513-8388

WEB: www.alfascientific.com • EMAIL: info@alfascientific.com