

Simple, Accurate Testing For Early Detection of Pregnancy



- FDA 510(k) Cleared for Professional Use
- CE Marked for Sale in Europe
- OTC Cleared
- CLIA Waived

INSTANT-VIEW[®] hCG Test

One-Step Pregnancy Immunoassay

For the detection of human chorionic gonadatropin (hCG) in human urine or serum as an aid in the early detection of pregnancy.



WHERE VALUE MEETS QUALITY IN DIAGNOSTICS™

INSTANT-VIEW[®] Simple, Accurate Testing For Early Detection Of Pregnancy



- Simple One-Step Test
- Over 99% Accurate
- Fast Results in 2-5 Minutes
- Convenient Room Temperature Storage
- Can Detect Pregnancy as Early as the First Day of a Missed Period
- Multiple Test Formats

Product Description	Kit Size	Detection Cut-Off	FDA Cleared	Waived	CE Mark
hCG Urine Cassette	25	25 mIU/mL	✓	✓	✓
hCG Urine Dip Strip	50	25 mIU/mL	✓	✓	✓
hCG Urine/Serum Cassette	25	25 mIU/mL	✓	Urine	✓
hCG Urine/Serum Dip Strip	50	25 mIU/mL	✓	Urine	✓
hCG Urine Midstream	1	25 mIU/mL	OTC	✓	Self-Testing

To order or for further information please contact one of our representatives today:



13200 Gregg Street, Poway, CA 92064

PH: 858-513-3888 • TOLL FREE: 877-204-5071 • FAX: 858-513-8388

WEB: www.alfascientific.com • EMAIL: info@alfascientific.com