

# Simple, Fast Testing To Help Predict The Best Time To Conceive



- FDA 510(k)  
Cleared for  
Professional or  
FDA 510(k) Exempt  
for OTC Use
- CLIA Waived
- CE Marked for Sale  
in Europe



*INSTANT-VIEW*<sup>®</sup> LH Test

One-Step Ovulation Immunoassay

For the detection of human Luteinizing Hormone (hLH) in human urine to predict ovulation.



*WHERE VALUE MEETS QUALITY IN DIAGNOSTICS™*

# INSTANT-VIEW<sup>®</sup> Simple, Accurate Testing To Help Predict The Best Time To Conceive



- Simple One-Step Test
- Over 99% Accurate
- Fast Results in 4-7 Minutes
- Convenient Room Temperature Storage
- Can Predict LH Surges 14 Days Before Period Starts
- Multiple Test Formats
- Made in the USA

Product Description	Kit Size	Cut-Off Level	FDA Cleared	Waived	CE Mark
LH Urine Cassette	25	20 mIU/mL	✓	✓	✓
LH Urine Dip Strip	50	20 mIU/mL	✓	✓	✓
LH Urine Midstream (OTC Use)	1	20 mIU/mL	Exempt		

*To order or for further information please contact one of our representatives today:*



13200 Gregg Street, Poway, CA 92064

PH: 858-513-3888 • TOLL FREE: 877-204-5071 • FAX: 858-513-8388

WEB: [www.alfascientific.com](http://www.alfascientific.com) • EMAIL: [info@alfascientific.com](mailto:info@alfascientific.com)