

PRESS RELEASE

Revised 11-12-2010

Alfa Officially Launches New CE Marked Oral-View® Test for Sale in Europe

Over the last several years the interest in Saliva testing, otherwise known as Oral Fluid testing, has grown dramatically. With its inception first in lab based testing in 2000, many rapid diagnostic test manufacturers saw an opportunity to develop tests that could provide results in as little as 5 minutes, much the same way urine-tests do. But none to date have been like this. With its unique, modern appearance and easy to use one-step design, Alfa's Oral-View® Drug of Abuse test is sure to be the next best thing in Oral Fluid screening.

The Oral-View® Saliva Multi-Drug of Abuse Test is a rapid, lateral flow diagnostic test designed to simultaneously detect up to 8 drugs of abuse in human saliva. This true one-step oral fluid testing device can be used to screen a donor and then itself be sent off to a laboratory for LC/MS confirmation. The Oral-View® test can screen for:

Abbreviation	Description	Cut-Off Level
AMP	Amphetamine	50 ng/mL
BZD	Benzodiazepine	20 ng/mL
COC	Cocaine	20 ng/mL
MET	Methamphetamine with Ecstasy (MDMA)	50 ng/mL
MOR/OPI	Morphine/Opiates	40 ng/mL
MTD	Methadone	30 ng/mL
PCP	Phencyclidine	10 ng/mL
THC	Marijuana	12 ng/mL

In many ways Oral Fluid (Saliva) testing has advantages over the traditional urine-based testing:

- Since collections can be observed, it is difficult for the donor to adulterate the specimen.
- Eliminates shy bladder issues that delay sample collections.
- It is a good indicator of recent drug usage. With urine testing, you must wait for the drug to be metabolized by the liver before it will show a preliminary positive result.
- Collections can be observed without having to worry about gender or donor privacy issues typically associated with urine testing.



13200 Gregg Street, Poway, CA 92064
PH: 858-513-3888 FX: 858-513-8388
TOLL FREE: 877-204-5071
Contact: Alison Adduono (aadduono@alfascientific.com)

PRESS RELEASE

- Cost effective as it eliminates the need for special collection facilities. Oral Fluid (Saliva) tests can be performed anytime - anywhere.

Important Notes: This product is available for export and is now CE Marked for sale in Europe. In the United States this product is for Forensic Use Only and this product is not for At Home Use or Over-the-Counter Sale. This product is not for use in any medical diagnosis or treatment in the United States.

For more information on our Oral-View® test or any of our other diagnostic tests please contact your Alfa sales or customer service representative today or visit our website at www.alfascientific.com.

About Alfa Scientific Designs, Inc.

Alfa Scientific Designs, Inc. is located in Poway, CA, and is a manufacturer and distributor of US made rapid diagnostic tests for Drugs of Abuse, Infectious Disease, Cardiac Markers and Fertility. Alfa offers a variety of test formats as well as OEM custom manufacturing capabilities.

Alfa-- "Where Value Meets Quality In Diagnostics."



13200 Gregg Street, Poway, CA 92064
PH: 858-513-3888 FX: 858-513-8388
TOLL FREE: 877-204-5071
Contact: Alison Adduono (aadduono@alfascientific.com)