

DRUGS OF ABUSE, INTEGRATED CUPS (URINE)

Additional Cup Configurations available – Contact us for details

PRODUCT	PRODUCT FORMAT	US FDA CLEARED	US CLIA WAIVED	CE MARKED	KIT CAT. NO	KIT SIZE	STOCK ITEM
AMP, COC, MET, OPI, THC	5 Drug Vision™ Cup	•		•	03-5736	25	•
AMP, COC, OPI, PCP, THC	5 Drug Vision™ Cup	•		•	03-5744	25	•
AMP, BZD, COC, MET, OPI, THC	6 Drug Vision™ Cup	•		•	03-5745	25	•
AMP, BAR, BZD, COC, MAMP, MDMA, MTD, OPI, OXY, PCP, PPX, TCA, THC	13 Drug Vision™ Cup	•		•	03-5752	25	•
AMP, BAR, BZD, COC, MET, MOR/OPI, MTD, PCP, TCA, THC, XTC	11 Drug Instant-View® SplitScreen Cup	•	•	•	03-4570	25	•

URINE COLLECTION CUPS

PRODUCT	SIZE	KIT CAT. NO	KIT SIZE	STOCK ITEM
Specimen Collection Cup with Flip Top	60 ml	55-3003	Bulk	•
Specimen Collection Cup with Screw Top Lid	60 ml	55-3004	Bulk	•
Optional Lid Insert	n/a	55-3006	Bulk	•

DRUGS OF ABUSE, MULTI DRUG PANEL, SALIVA - ORAL-VIEW® BRAND

First Generation - One Strip Design

PRODUCT	PRODUCT FORMAT	US FDA CLEARED	US CLIA WAIVED	CE MARKED	KIT CAT. NO	KIT SIZE	STOCK ITEM
THC (For Export)	1 Drug Cassette			•	03-4203	25	
COC, OPI, THC (For Export)	3 Drug Cassette			•	03-4205	25	
AMP, BZD, COC, OPI, THC (For Export)	5 Drug Cassette			•	03-4210	25	•
AMP, BZD, COC, OPI, THC (US - "For Forensic or Investigational Use Only")	5 Drug Cassette				03-4231	25	•

DRUGS OF ABUSE, MULTI DRUG PANEL, SALIVA - ORAL-VIEW® BRAND

Second Generation - Two Strip Design

PRODUCT	PRODUCT FORMAT	US FDA CLEARED	US CLIA WAIVED	CE MARKED	KIT CAT. NO	KIT SIZE	STOCK ITEM
AMP, COC, OPI, THC, MET/XTC (For Export)	5 Drug Cassette			•	03-4261	25	•
AMP, COC, OPI, THC, MET/XTC (US - "For Forensic or Investigational Use Only")	5 Drug Cassette				03-4262	25	
AMP, BZD, COC, OPI, THC, MET/XTC (For Export)	6 Drug Cassette			•	03-4271	25	•
AMP, BZD, COC, OPI, THC, MET/XTC (US - "For Forensic or Investigational Use Only")	6 Drug Cassette				03-4272	25	•
AMP, BZD, COC, OPI, PCP, MET/XTC, MTD, THC (For Export)	8 Drug Cassette			•	03-4281	25	
AMP, BZD, COC, OPI, PCP, MET/XTC, MTD, THC (US - "For Forensic or Investigational Use Only")	8 Drug Cassette				03-4282	25	•

DRUGS OF ABUSE, MISCELLANEOUS TESTS (URINE)

Made in China

PRODUCT	PRODUCT FORMAT	DESCRIPTION	KIT CAT. NO	KIT SIZE	STOCK ITEM
Adulteration Strips	Dip Strip	Creatinine, Nitrite, Glutaraldehyde, pH, Specific Gravity, PCC/Oxidants	03-0012	25	•
Cotinine	Cassette	Cotinine (Nicotine) - <i>For Forensic Use Only. Not for sale in Europe.</i>	03-0050	40	

DRUGS OF ABUSE, SINGLE TESTS, INSTANT-VIEW® BRAND (URINE)

PRODUCT	PRODUCT FORMAT	SENSITIVITY	DESCRIPTION	US FDA CLEARED	US CLIA WAIVED	CE MARKED	KIT CAT. NO	KIT SIZE	STOCK ITEM
AMP	Cassette	1000 ng/ml	Amphetamine	•		•	03-3012	25	•
	Dip Strip			•		•	03-3017	50	•
	Dip Cassette			•	•		03-3012QW	25	•
AMP300	Cassette	300 ng/ml	Amphetamine	•		•	03-3014	25	
	Dip Strip			•		•	03-3019	50	
AMP500	Cassette	500 ng/ml	Amphetamine			•	03-3010	25	
	Dip Strip					•	03-3018	50	
BAR	Cassette	200 ng/ml	Barbiturate	•		•	03-3152	25	•
	Dip Strip			•		•	03-3157	50	•
	Dip Cassette			•	•		03-3152QW	25	•
BZD	Cassette	300 ng/ml	Benzodiazepine	•		•	03-3132	25	•
	Dip Strip			•		•	03-3137	50	•
	Dip Cassette			•	•		03-3132QW	25	•
BZD200	Cassette	200 ng/ml	Benzodiazepine				03-3143	25	
	Dip Strip						n/a	50	
BUP	Cassette	10 ng/ml	Buprenorphine	•		•	03-2842	25	•
	Dip Strip			•		•	03-2847	50	
COC	Cassette	300 ng/ml	Cocaine	•		•	03-3052	25	•
	Dip Strip			•		•	03-3057	50	•
	Dip Cassette			•	•		03-3052QW	25	•
COC150	Cassette	150 ng/ml	Cocaine	•		•	03-3054	25	
	Dip Strip			•		•	03-3058	50	
EDDP	Cassette	100 ng/ml	EDDP (Methadone Metabolite)				03-4842	25	
	Dip Strip						03-4847	50	
KET	Cassette	1000 ng/ml	Ketamine				03-4802	25	
	Dip Strip						03-4807	50	
METH	Cassette	1000 ng/ml	Methamphetamine	•		•	03-3032	25	•
	Dip Strip			•		•	03-3037	50	•
	Dip Cassette			•	•		03-3032QW	25	•
METH500	Cassette	500 ng/ml	Methamphetamine	•		•	03-3033	25	
	Dip Strip			•		•	03-3038	50	
METH300	Cassette	300 ng/ml	Methamphetamine	•		•	03-3030	25	
	Dip Strip			•		•	03-3035	50	
MOR/OPI	Cassette	2000 ng/ml	Morphine/Opiates	•		•	03-3092	25	•
	Dip Strip			•		•	03-3097	50	•
	Dip Cassette			•	•		03-3092QW	25	•
MOR/OPI 300	Cassette	300 ng/ml	Morphine/Opiates	•		•	03-3093	25	•
	Dip Strip			•		•	03-3098	50	•
	Dip Cassette			•	•		03-3093QW	25	•
MTD	Cassette	300 ng/ml	Methadone	•		•	03-3172	25	•
	Dip Strip			•		•	03-3177	50	•
	Dip Cassette			•	•		03-3172QW	25	•
OXY100	Cassette	100 ng/ml	Oxycodone	•		•	03-2862	25	
	Dip Strip			•		•	03-2867	50	
	Dip Cassette			•	•		03-2862QW	25	•
OXY300	Cassette	300 ng/ml	Oxycodone	•			03-2863	25	
	Dip Strip			•			03-2866	50	
	Dip Cassette			•	•		03-2863QW	25	
PCP	Cassette	25 ng/ml	Phencyclidine	•		•	03-3112	25	•
	Dip Strip			•		•	03-3117	50	•
	Dip Cassette			•	•		03-3112QW	25	•
PPX	Cassette	300 ng/ml	Propoxyphene	•		•	03-2932	25	•
	Dip Strip			•		•	03-2937	50	
TCA	Cassette	1000 ng/ml	Tricyclic Antidepressants	•		•	03-2972	25	•
	Dip Strip			•		•	03-2977	50	
	Dip Cassette			•	•		03-2972QW	25	•
THC	Cassette	50 ng/ml	Marijuana	•		•	03-3072	25	•
	Dip Strip			•		•	03-3077	50	•
	Dip Cassette			•	•		03-3072QW	25	•
TRA	Cassette	200 ng/ml	Tramadol			•	03-4862	25	
	Dip Strip					•	03-4867	50	
XTC	Cassette	500 ng/ml	MDMA (Ecstasy)	•		•	03-2952	25	•
	Dip Strip			•		•	03-2957	50	
	Dip Cassette			•	•		03-2952QW	25	•

DRUGS OF ABUSE, MULTI DRUG PANEL CONFIGURATIONS - INSTANT-VIEW® BRAND (URINE)

Additional Multi Panel Configurations available – Contact us for details

Urine Adulteration (Creatinine, Nitrite, Glutaraldehyde, pH, Specific Gravity, Oxidants) are specimen adulteration controls

PRODUCT	PRODUCT FORMAT	US FDA CLEARED	US CLIA WAIVED	CE MARKED	KIT CAT. NO	KIT SIZE	STOCK ITEM
COC, THC	2 Drug Dual Format Cassette	•		•	03-3270	25	•
COC, MET, THC	3 Drug Dual Format Cassette	•		•	03-3371	25	•
COC, MOR/OPI, THC	3 Drug Dual Format Cassette	•		•	03-3373	25	•
AMP, COC, THC	3 Drug Dual Format Cassette	•		•	03-3375	25	•
COC, MOR/OPI 300, THC	3 Drug Dual Format Cassette	•		•	03-3377	25	
AMP, COC, MOR/OPI 300, THC w/ADULTERATION	4 Drug Dual Format Cassette	•		•	03-3406	25	
AMP, COC, MOR/OPI, THC	4 Drug Dual Format Cassette	•		•	03-3458	25	•
COC, MET, MOR/OPI 300, THC	4 Drug Dual Format Cassette	•		•	03-3462	25	•
COC, MET, MOR/OPI, THC	4 Drug Dual Format Cassette	•		•	03-3474	25	
BZD, COC, MOR/OPI, MTD	4 Drug CLIA Dip Card	•	•	•	03-3427QW	25	
COC, MET, MOR/OPI, PCP, THC w/ADULTERATION	5 Drug Dual Format Cassette	•		•	03-3502	25	•
COC, MET, MOR/OPI, THC, XTC	5 Drug CLIA Dip Card	•	•	•	03-3531QW	25	•
BZD, COC, MET, MOR/OPI, THC	5 Drug Dual Format Cassette	•		•	03-3571	25	•
BZD, COC, MET, MOR/OPI, THC	5 Drug CLIA Dip Card	•	•	•	03-3571QW	25	•
AMP, COC, MOR/OPI, PCP, THC	5 Drug Dual Format Cassette	•		•	03-3576	25	•
AMP, COC, MOR/OPI, PCP, THC	5 Drug CLIA Dip Card	•	•	•	03-3576QW	25	•
AMP, COC, MET, MOR/OPI, THC	5 Drug Dual Format Cassette	•		•	03-3577	25	•
AMP, COC, MET, MOR/OPI, THC	5 Drug CLIA Dip Card	•	•	•	03-3577QW	25	•
COC, MET, MOR/OPI, PCP, THC	5 Drug Dual Format Cassette	•		•	03-3578	25	•
COC, MET, MOR/OPI, PCP, THC	5 Drug CLIA Dip Card	•	•	•	03-3578QW	25	•
COC, MET, MOR/OPI 300, PCP, THC	5 Drug Dual Format Cassette	•		•	03-3583	25	•
AMP, BZD, COC, MET, MOR/OPI, THC	6 Drug Dual Format Cassette	•		•	03-3635	25	•
AMP, BZD, COC, MET, MOR/OPI, THC	6 Drug CLIA Dip Card	•	•	•	03-3635QW	25	•
AMP, BAR, BZD, COC, MOR/OPI 300, THC	6 Drug Dual Format Cassette	•		•	03-3669	25	
AMP, COC, MET, MOR/OPI, PCP, THC	6 Drug Dual Format Cassette	•		•	03-3679	25	•
AMP, COC, MET, MOR/OPI, PCP, THC	6 Drug CLIA Dip Card	•	•	•	03-3679QW	25	•
AMP, BAR, BZD, COC, MET, MOR/OPI, THC	7 Drug Dual Format Cassette	•		•	03-3761	20	
AMP, BAR, BZD, COC, MET, MOR/OPI, THC	7 Drug CLIA Dip Card	•	•	•	03-3761QW	20	
AMP, BAR, BZD, COC, MET, MOR/OPI, MTD, PCP, THC	9 Drug Dual Format Cassette	•		•	03-3997	20	•
AMP, BAR, BZD, COC, MET, MOR/OPI, MTD, PCP, THC	9 Drug CLIA Dip Card	•	•	•	03-3997QW	20	•
AMP, BAR, BZD, COC, MET, MOR/OPI, MTD, PCP, TCA, THC	10 Drug Dual Format Cassette	•		•	03-3910	20	•
AMP, BAR, BZD, COC, MET, MOR/OPI, MTD, PCP, TCA, THC	10 Drug CLIA Dip Card	•	•	•	03-3910QW	20	•
AMP, BAR, BZD, COC, MET, MOR/OPI, MTD, PCP, THC, XTC	10 Drug CLIA Dip Card	•	•	•	03-3914QW	20	
AMP, BAR, BZD, COC, MET, MOR/OPI 300, MTD, PCP, THC, XTC	10 Drug Dual Format Cassette	•		•	03-3917	20	•
AMP, BAR, BZD, COC, MET, MOR/OPI, MTD, OXY, PCP, THC	10 Drug Dual Format Cassette	•		•	03-3929	20	
AMP, BAR, BZD, COC, MET, MOR/OPI, OXY, PPX, THC, XTC	10 Drug Dual Format Cassette	•		•	03-3922	20	•
AMP, BAR, BZD, COC, MET, MOR/OPI, MTD, PCP, TCA, THC, XTC	11 Drug CLIA Dip Card	•	•	•	03-3940QW	20	•
AMP, BAR, BZD, COC, MET, MOR/OPI, MTD, PCP, PPX, TCA, THC, XTC	12 Drug Dual Format Cassette	•		•	03-3951	20	•
AMP, BAR, BZD, COC, MET, MOR/OPI 300, MTD, PCP, PPX, TCA, THC, XTC	12 Drug Dual Format Cassette	•		•	03-3953	20	•
AMP, BAR, BZD, COC, MET, MOR/OPI 300, MTD, OXY, PCP, TCA, THC, XTC	12 Drug CLIA Dip Card	•	•	•	03-3965QW	20	•

FERTILITY - INSTANT-VIEW® BRAND

PRODUCT	SPECIMEN	PRODUCT FORMAT	SENSITIVITY	DESCRIPTION	US FDA CLEARED	US CLIA WAIVED	CE MARKED	KIT CAT. NO	KIT SIZE	STOCK ITEM
hCG	Urine	Cassette	25 mIU	Human Chorionic Gonadotropin (Pregnancy)	•	•	•	02-2482	25	•
	Urine	Dip Strip			•	•	•	02-2487	50	•
	Urine	Midstream			OTC	•	Self-testing	07-2480	1	Bulk
	Urine/Serum	Cassette			•	Urine	•	02-2483	25	•
	Urine/Serum	Dip Strip			•	Urine	•	02-2488	50	•
hCG II	Urine/Serum	Cassette	10 mIU	Pregnancy			•	02-2481	25	•
hCG Ref.	Urine/Serum	Cassette	10 & 20 mIU	Pregnancy			•	02-2485	25	
LH	Urine	Cassette	20 mIU	Luteinizing Hormone (Ovulation)	•	•	•	02-2542	25	•
	Urine	Dip Strip			•	•	•	02-2547	50	•
	Urine	Midstream			OTC Exempt	•		07-2540	1	Bulk
FSH	Urine	Cassette	40 mIU	Follicle Stimulating Hormone (Menopause)	Exempt	•	•	02-2562	25	
	Urine	Dip Strip			Exempt	•	•	02-2567	50	
	Urine	Midstream			Exempt	•		07-2560	1	

FERTILITY - GENERIC

Assembled in China

PRODUCT	SPECIMEN	PRODUCT FORMAT	SENSITIVITY	DESCRIPTION	US FDA CLEARED	US CLIA WAIVED	CE MARKED	KIT CAT. NO	KIT SIZE	STOCK ITEM
hCG	Urine	Cassette	25 mIU	Human Chorionic Gonadotropin (Pregnancy)	•	•		07-2492A	Bulk	•
	Urine	Dip Strip			•	•		07-2497A	Bulk	•
	Urine/Serum	Cassette			•	Urine		07-2493A	Bulk	•

INFECTIOUS DISEASES - INSTANT-VIEW® BRAND

PRODUCT	SPECIMEN	PRODUCT FORMAT	DESCRIPTION	US FDA CLEARED	US CLIA WAIVED	CE MARKED	KIT CAT. NO	KIT SIZE	STOCK ITEM
H. pylori	Serum	Cassette	Helicobacter pylori	•		•	05-9473	25	•
	Serum	Dip Strip		•		•	05-9478	50	
	Serum/Whole Blood	Cassette		•		•	05-9494	25	•
	Whole Blood	Cassette		•	•	•	05-9475	25	•
HBsAg	Serum	Cassette	Hepatitis B Virus Surface Antigen				05-9843	25	•
	Serum	Dip Strip					05-9848	50	•
HCV	Serum	Cassette	Hepatitis C Virus Antibody				05-9853	25	•
Strep A	Throat Swab	Cassette	Group A Streptococci			•	05-9301	25	
	Throat Swab	Dip Strip				•	05-9302	50	

CANCER MARKERS - INSTANT-VIEW® BRAND

PRODUCT	SPECIMEN	PRODUCT FORMAT	SENSITIVITY	DESCRIPTION	US FDA CLEARED	US CLIA WAIVED	CE MARKED	KIT CAT. NO	KIT SIZE	STOCK ITEM
PSA Semi-quantitative	Serum	Cassette	1 & 4 ng/ml	Prostate Specific Antigen				05-7772	25	•
	Serum/Whole Blood	Cassette						05-7771	25	•

CARDIAC MARKERS - INSTANT-VIEW® BRAND

PRODUCT	SPECIMEN	PRODUCT FORMAT	SENSITIVITY	DESCRIPTION	US FDA CLEARED	US CLIA WAIVED	CE MARKED	KIT CAT. NO	KIT SIZE	STOCK ITEM
TNI	Serum	Cassette	1.5 ng/ml free or 5 ng/ml complex cTnl	Cardiac Troponin I	•		•	05-8572	25	•
	Serum/Whole Blood	Cassette			•		•	05-8575	25	•
3 in 1	Serum/Whole Blood	Cassette	See package insert	Troponin I CK-MB Myoglobin			•	05-8503	25	

Products without US FDA clearance cannot be sold in USA. Forensic Use Only products may also be sold in US.
Products without CE Mark cannot be sold in EU.