

# Simple, Fast Testing Find Out If They Have Been Drinking



**NEW!**

## Oral-View<sup>®</sup> Saliva One-Step Alcohol Assay

For the detection of alcohol easily and quickly.

- 510K Cleared and CLIA-Waived
- Accurate and Reliable
- Convenient
- Results in 2 Minutes



**WHERE VALUE MEETS QUALITY IN DIAGNOSTICS™**

# Oral-View® Simple, Accurate Testing For Detecting Alcohol in Saliva



- Quick and Easy
- Semi-Quantitative Results
- Suitable for use in homes, schools, workplace, military, law enforcement and other government agencies
- Non-invasive specimen collection
- No equipment needed
- Cost efficient and effective

Features	Details
Product Name	Saliva Alcohol Test Strips
Catalog Number	03-0007
Regulatory Approval	510K Cleared & CLIA-Waived
Technology	Enzymatic Reaction
Specimen	Saliva
Blood Alcohol Concentration (BAC)	0.02% to 0.30%
Measuring Units	% BAC
Reading Time	2 minutes
Reading Stability	3 minutes
Storage Temperature	2-27°C (36-81°F)
Shelf Life	1 year
Pouch Dimension (L X W)	13.5 x 4.7 cm (5.3" x 1.9")
Kit Box Dimension (L X W X H)	18.0 x 10.5 x 5.0 cm (7.1" x 4.1" x 2.0")
Package Size (# of Tests)	25

To order or for further information please contact one of our representatives today:



13200 Gregg Street, Poway, CA 92064

PH: 858-513-3888 • TOLL FREE: 877-204-5071 • FAX: 858-513-8388

WEB: [www.alfascientific.com](http://www.alfascientific.com) • EMAIL: [info@alfascientific.com](mailto:info@alfascientific.com)