


Simple, Fast Testing For Detecting Drugs Of Abuse in Saliva



NEW!

- CE Marked for Sale in Europe
- For Forensic Use Only in US Market



Patented
US Patent No. D626249

Oral-View[®] Saliva One-Step Drugs of Abuse Immunoassay

For the detection of up to 8 illegal or prescription drugs in human saliva to determine illicit use.



“WHERE VALUE MEETS QUALITY IN DIAGNOSTICS”™

Oral-View[®] Simple, Accurate Testing For Detecting Drugs Of Abuse in Saliva



- Results in About 5 Minutes
- One Step Testing
- Easy Collection Method
- Minimal Sample Required
- Device Can Be Sent For Positive Result Confirmation
- Also Available in Single Drug Formats
- Made in USA

DRUG ANALYTE CHOICES		
Product	Description	Cut Off
AMP	Amphetamine	50 ng/mL
BZD	Benzodiazepine	20 ng/mL
COC	Cocaine	20 ng/mL
MET	Methamphetamine & Ecstasy (MDMA)	50 ng/mL
OPI	Morphine/Opiate	40 ng/mL
MTD	Methadone	30 ng/mL
PCP	Phencyclidine	10 ng/mL
THC	Marijuana	12 ng/mL

To order or for further information please contact one of our representatives today:



13200 Gregg Street, Poway, CA 92064

PH: 858-513-3888 • TOLL FREE: 877-204-5071 • FAX: 858-513-8388

WEB: www.alfascientific.com • EMAIL: info@alfascientific.com

00-0125 F

“WHERE VALUE MEETS QUALITY IN DIAGNOSTICS”™