

Simple, Accurate Testing We Have The *VISION*TM - Do You?

NEW!

- FDA 510(k)
Cleared for
Professional Use
- CE Marked for Sale
in Europe



*VISION*TM DOA Test Cup

One-Step Drugs of Abuse Immunoassay

For the detection of illegal or prescription drugs in human urine to drug abuse.
Ideal for pre-employment and workplace testing.




WHERE VALUE MEETS QUALITY IN DIAGNOSTICSTM

VISION™ Simple, Accurate DOA Testing

DRUG ANALYTE CHOICES				
Product	Description	Cut Off	FDA Cleared	CE Mark
AMP	Amphetamine	1000 ng/mL	✓	✓
AMP300	Amphetamine	300 ng/mL	✓	✓
BAR	Barbiturate	300 ng/mL	✓	✓
BZD	Benzodiazepine	300 ng/mL	✓	✓
BUP	Buprenorphine	10 ng/mL	✓	✓
COC	Cocaine	300 ng/mL	✓	✓
COC150	Cocaine	150 ng/mL	✓	✓
MET	Methamphetamine	1000 ng/mL	✓	✓
MET500	Methamphetamine	500 ng/mL	✓	✓
MDMA	Ecstasy/MDMA	500 ng/mL	✓	✓
MTD	Methadone	300 ng/mL	✓	✓
OPI	Morphine/Opiate	2000 ng/mL	✓	✓
MOP	Morphine/Opiate	300 ng/mL	✓	✓
OXY	Oxycodone	100 ng/mL	✓	✓
PCP	Phencyclidine	25 ng/mL	✓	✓
PPX	Propoxyphene	300 ng/mL	✓	✓
TCA	Tricyclic Antidepressant	1000 ng/mL	✓	✓
THC	Marijuana	50 ng/mL	✓	✓

- One Step Operation
- Results in 5 Minutes or Less
- Results Stable for 60 Minutes
- Sample Access Port for Easy Confirmation
- Leakproof Cap Locking Feature
- Stackable for Easy, Less Costly Shipping



To order or for further information please contact one of our representatives today:



13200 Gregg Street, Poway, CA 92064

PH: 858-513-3888 • TOLL FREE: 877-204-5071 • FAX: 858-513-8388

WEB: www.alfascientific.com • EMAIL: info@alfascientific.com