



Instant-view ^{plus}™ Simple Cup



**RESULTS AS
EARLY AS
2 MINUTES**

- One Step Collection and Test
- Patent Protected
- Maximum Configuration
13 Drug Analytes

Classification	# of Drug Analytes
Professional	20
CLIA	13
CE	25
Forensic	27



Patent #9,921,217

Accuracy at the Speed of ALFA




Instant-view ^{plus}™

AVAILABLE ANALYTES

Abbreviation	Analytes	Cutoff (ng/mL)	PRO	CLIA	CE	Forensic Use
AMP	Amphetamine	1000	√	√	√	√
AMP 300	Amphetamine	300	√		√	√
AMP 500	Amphetamine	500			√	√
BAR	Barbiturates	200	√	√	√	√
BUP/NBUP	Buprenorphine/Norbuprenorphine	10	√	√	√	√
BZD	Benzodiazepine	200	√	√	√	√
COC	Cocaine	300	√	√	√	√
COC 150	Cocaine	150	√		√	√
EtG	Ethyl Glucuronide	300			√	√
FTL	Fentanyl	20			√	√
FTL	Fentanyl	200				√
HRN	6-Acetylmorphine	10				√
K2	Synthetic Cannabis	50			√	√
MET	Methamphetamine	1000	√	√	√	√
MET 300	Methamphetamine	300	√		√	√
MET 500	Methamphetamine	500	√		√	√
MOR	Morphine/Opiates	2000	√	√	√	√
MOR 300	Morphine/Opiates	300	√		√	√
MTD	Methadone	300	√	√	√	√
OXY	Oxycodone	100	√		√	√
OXY 300	Oxycodone	300	√	√	√	√
PCP	Phencyclidine	25	√	√	√	√
PPX	Propoxyphene	300	√		√	√
TCA	Tricyclics	1000	√	√	√	√
THC	Marijuana/Hashish	50	√	√	√	√
TRA	Tramadol	200			√	√
XTC	MDMA or Ecstasy	500	√	√	√	√



Read results at five (5) minutes.
Results stable up to sixty (60) minutes.



SAMPLE ONLY - Color chart NOT Exact