



Instant-View™
rapid medical tests

DRUG TESTS



US Patents
7.052.651 B2
D559 384 S
D529 170 S

CE Marked

US FDA **510(k)** OTC Cleared, **CLIA** Waived

- Test Up to 14 Drugs at Once
- Accuracy Up to 99.8%
- State-of-the-Art Design

Speed.
Accuracy.
Excellence.

Developed & Manufactured in the USA

US FDA Regulated • ISO Certified



Instant-View™
rapid medical tests

Specimen

✓ Urine

Drug Tests

Up to 14 Analytes in One Test • Adulteration Controls Available

DRUG ANALYTE CHOICES

Product	Description	Sensitivity (ng/ml)
AMP	Amphetamine	1000
AMP300	Amphetamine	300
BAR	Barbiturate	200
BZD	Benzodiazepine	300
BUP	Buprenorphine	10
COC	Cocaine	300
COC150	Cocaine	150
MET	Methamphetamine	1000
MET300	Methamphetamine	300
MET500	Methamphetamine	500
MOR	Morphine / Opiate Metabolite	2000
MOR/300	Morphine / Opiate Metabolite	300
MTD	Methadone	300
OXY	Oxycodone	300
OXY100	Oxycodone	100
PCP	Phencyclidine	25
PPX	Propoxyphene	300
TCA	Tricyclic Antidepressants	1000
THC	Marijuana	50
XTC	Ecstasy / MDMA	500

FORENSIC & RESEARCH USE ONLY - IN THE USA

EDDP	Methadone Metabolite	100
FTL	Fentanyl	200
K2	Synthetic Marijuana	50
KET	Ketamine	1000
TRA	Tramadol	200

Product Formats

- Dip Stick, Dip Card, Cassette, Cup

Product Registration Status:

- US FDA Professional / OTC **510(k)** Cleared
- **CLIA** Waived
- **CE** Marked
- Forensic / Research

Made in
the USA

To Order • For Further Information • For Live Tech Support

TOLL FREE: 877-204-5071

www.alfascientific.com • info@alfascientific.com



Rapid Medical Tests

Alfa Scientific Designs, Inc.
13200 Gregg Street
Poway, CA 92064 – USA