



## HEALTH TESTS

- Infectious Disease

## H. (HELICOBACTOR) PYLORI TEST

- Over 95.1% Sensitive
- Over 92.8% Specific
- Fast Results in 4-7 Minutes
- Room Temperature Storage
- Specific for Toxic Strains

Speed.
Accuracy.
Excellence.

**Developed & Manufactured in the USA** 

**US FDA Regulated • ISO Certified** 





## **Health Tests**

Gastrointestinal Disorder

HEALTH TESTS - INFECTIOUS DISEASE							
Specimen: Serum / Whole Blood							
Test Format: Cassette							
PRODUCT	DESCRIPTION	FORMAT	SPECIMEN	SENSITIVITY ng/ml	REGISTRATION STATUS		
					510(k) Professional	Forensic/ Research	CE Marked
			Serum	- Taria Charles	<b>~</b>		<b>✓</b>
H. pylori	Helicobacter pylori	Cassette	Serum Whole Blood	Toxic Strains	~		•

## Intended use - H. pylori

The H. pylori test is a rapid lateral flow, qualitative immunoassay. It is intended to detect the presence of antibodies against Helicobacter pylori in human blood or serum.



To Order • For Further Information • For Live Tech Support

TOLL FREE: 877-204-5071

www.alfascientific.com • info@alfascientific.com



Alfa Scientific Designs, Inc. 13200 Gregg Street Poway, CA 92064 – USA